
GlaucoScan

Overview

This project is a deep learning-based desktop application designed to assist in the preliminary inspection of **Fundus Camera Images** for signs of **Glaucoma**.

It utilizes a three-stage computer vision pipeline:

1. **Preprocessing:** This included background removal to eliminate black borders, contrast enhancement using gamma correction and other image augmentation.
2. **UNet Segmentation:** Accurately localizes and segments the optic disc and cup region.
3. **EfficientNet-B3 Classification:** Classifies the cropped region to distinguish between **Non-Glaucoma** and potential **Glaucoma** cases.

The system is built using **PyTorch** and deployed with a user-friendly **PyQt5 GUI**, allowing researchers and medical professionals to perform local, real-time image analysis.



Features

- **Two-Stage Preprocessing:** Implements robust preprocessing and UNet segmentation to ensure the classification model focuses only on the critical optic disc region, maximizing diagnostic accuracy.
- **High-Performance Backbone:** Utilizes a fine-tuned **EfficientNet-B3** network for fast and reliable feature extraction and classification.
- **Fundus Image Optimized:** The segmentation and preprocessing steps are specifically tailored for varying qualities and sizes of fundus camera images, effectively removing dark backgrounds and centering on the region of interest.
- **Confidence Metrics:** Displays precise probability scores for both the "Glaucoma" and "Non-Glaucoma" categories.
- **Clean Desktop UI:** A fully featured, self-contained desktop interface built with PyQt5, optimizing ease of use for local image analysis.



Installation and Setup

Here we provide two ways: python scripts and software

For python scripts

Prerequisites

- Python 3.8+
- NVIDIA GPU (Recommended for faster inference, otherwise runs on CPU)

Environment Setup

1. Clone the repository:

```
1. git clone https://github.com/Brokengm/BIA_group_work.git
2. cd final_version
```

2. Create and activate the environment (recommend conda, if you don't have conda, you can search installation guidance online):

```
1. # Using conda
2. conda create -n glaucoma python=3.9
3. conda activate glaucoma
```

3. Install dependencies:

```
1. pip install torch torchvision torchaudio --index-url
   https://download.pytorch.org/whl/cu118 # Use appropriate CUDA version
2. pip install -r requirements.txt
3. # The required packages include: PyQt5, numpy, Pillow, timm, scikit-image
```

4. Files you need:

For the application to run, you must download these files and place them directly into the same directory:

- `384_unet.pth` UNet Segmentation Model
- `efficientnet.pth` EfficientNet Classification Model
- `app_gui.py` Main application, GUI (PyQt5) logic, and prediction pipeline entry point.
- `model.py` Contains UNet and CombinedModel (EfficientNet) class definitions and model loading logic.
- `preprocess.py` Contains all image processing functions (cropping, CLAHE, segmentation, etc.).
- `welcome.jpeg` The picture of main page.

If your system is MacOS, you can try our software

Just download GlaucoScan.zip, and decompress.



Usage

Running the Application

1. Ensure you have completed the **Installation and Setup** steps.
2. Run the main application script:

```
1. python app_gui.py
```

Analysis Steps

1. **Welcome Screen:** Click **“Let’s get started”** .
2. **Load Image:** Click the **“Add image”** button to select a fundus image (, , etc.).
The image will appear in the left panel.
3. **Analyze:** Click **“Run analysis”** . The system will perform:
 - Preprocessing (border removal, contrast enhancement).
 - Optic Disc Segmentation (using UNet).
 - Optic Disc Cropping (centering on the segmented region).
 - Glaucoma Classification (using EfficientNet-B3).
4. **View Results:** The classification result, including the predicted class (or) and the precise probability scores, will be displayed in the right panel.



Disclaimer

This software is intended for course work and educational demonstration only. It must NOT be used for real medical diagnosis or self-assessment of glaucoma. Always consult a qualified medical professional for health concerns.